

Antibodies Product Datasheet

Product Name:	FLNB Antibody
Catalog Number:	C49421
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FLNB
Concentration:	1 mg/ml
Gene Name:	FLNB FLN1L FLN3 TABP TAP
Protein Name:	Filamin-B (FLN-B) (ABP-278) (ABP-280 homolog) (Actin-binding-like protein) (Beta-filamin) (Filamin homolog 1) (Fh1) (Filamin-3) (Thyroid autoantigen) (Truncated actin-binding protein) (Truncated ABP)
Other Name:	
Human Swiss Prot No.:	O75369
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FLNB at Human
Storage and Stability:	Store at -20° C/1 year