

Antibodies Product Datasheet

Product Name:	FEM1B Antibody
Catalog Number:	C49426
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FEM1B
Concentration:	1 mg/ml
Gene Name:	FEM1B F1AA KIAA0396
Protein Name:	Protein fem-1 homolog B (FEM1b) (FEM1-beta) (Fem-1-like death receptor-binding protein alpha) (Fem-1-like in apoptotic pathway protein alpha) (F1A-alpha)
Other Name:	
Human Swiss Prot No.:	Q9UK73
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FEM1B at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year