

Antibodies Product Datasheet

Product Name:	GTPBP9 Antibody
Catalog Number:	C49427
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GTPBP9
Concentration:	1 mg/ml
Gene Name:	OLA1 GTPBP9 PRO2455 PTD004
Protein Name:	Obg-like ATPase 1 (EC 3.6.3.-) (DNA damage-regulated overexpressed in cancer 45) (DOC45) (GTP-binding protein 9)
Other Name:	
Human Swiss Prot No.:	Q9NTK5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GTPBP9 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year