

Antibodies Product Datasheet

Product Name:	EXD2 Antibody
Catalog Number:	C49437
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human EXD2
Concentration:	1 mg/ml
Gene Name:	EXD2 C14orf114 EXDL2
Protein Name:	Exonuclease 3'-5' domain-containing protein 2 (Exonuclease 3'-5' domain-like-containing protein 2)
Other Name:	
Human Swiss Prot No.:	Q9NVH0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of EXD2 at Human
Storage and Stability:	Store at -20° C/1 year