

Antibodies Product Datasheet

Product Name:	HVCN1 Antibody
Catalog Number:	C49438
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HVCN1
Concentration:	1 mg/ml
Gene Name:	HVCN1 VSOP UNQ578/PRO1140
Protein Name:	Voltage-gated hydrogen channel 1 (Hydrogen voltage-gated channel 1) (HV1) (Voltage sensor domain-only protein)
Other Name:	
Human Swiss Prot No.:	Q96D96
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HVCN1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year