

Antibodies Product Datasheet

Product Name:	CNKS3 Antibody
Catalog Number:	C49451
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CNKS3
Concentration:	1 mg/ml
Gene Name:	CNKS3 MAGI1
Protein Name:	Connector enhancer of kinase suppressor of ras 3 (Connector enhancer of KSR 3) (CNK homolog protein 3) (CNK3) (CNKS3 family member 3) (Maguin-like protein)
Other Name:	
Human Swiss Prot No.:	Q6P9H4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CNKS3 at Human, Mouse, Rat
Storage and Stability:	Store at -20° C/1 year