

Antibodies Product Datasheet

Product Name:	BBS1 Antibody
Catalog Number:	C49457
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human BBS1
Concentration:	1 mg/ml
Gene Name:	BBS1 BBS2L2
Protein Name:	Bardet-Biedl syndrome 1 protein (BBS2-like protein 2)
Other Name:	
Human Swiss Prot No.:	Q8NFJ9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BBS1 at Human
Storage and Stability:	Store at -20° C/1 year