

Antibodies Product Datasheet

Product Name:	Harmonin Antibody
Catalog Number:	C49468
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human Harmonin
Concentration:	1 mg/ml
Gene Name:	USH1C AIE75
Protein Name:	Harmonin (Antigen NY-CO-38/NY-CO-37) (Autoimmune enteropathy-related antigen AIE-75) (Protein PDZ-73) (Renal carcinoma antigen NY-REN-3) (Usher syndrome type-1C protein)
Other Name:	
Human Swiss Prot No.:	Q9Y6N9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Harmonin at Human, Mouse
Storage and Stability:	Store at -20° C/1 year