

Antibodies Product Datasheet

Product Name:	LIN9 Antibody
Catalog Number:	C49477
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human LIN9
Concentration:	1 mg/ml
Gene Name:	LIN9 BARA TGS
Protein Name:	Protein lin-9 homolog (HuLin-9) (hLin-9) (Beta subunit-associated regulator of apoptosis) (TUDOR gene similar protein) (Type I interferon receptor beta chain-associated protein) (pRB-associated protein)
Other Name:	
Human Swiss Prot No.:	Q5TKA1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LIN9 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year