

Antibodies Product Datasheet

Product Name:	AKR7A3 Antibody
Catalog Number:	C49488
Reactivity:	Human, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human AKR7A3
Concentration:	1 mg/ml
Gene Name:	AKR7A3 AFAR2
Protein Name:	Aflatoxin B1 aldehyde reductase member 3 (EC 1.-.-.) (AFB1 aldehyde reductase 2) (AFB1-AR 2)
Other Name:	
Human Swiss Prot No.:	O95154
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of AKR7A3 at Human,Rat
Storage and Stability:	Store at -20° C/1 year