

Antibodies Product Datasheet

Product Name:	ARD1B Antibody
Catalog Number:	C49489
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ARD1B
Concentration:	1 mg/ml
Gene Name:	NAA11 ARD1B ARD2
Protein Name:	N-alpha-acetyltransferase 11 (EC 2.3.1.88) (N-terminal acetyltransferase complex ARD1 subunit homolog B) (hARD2) (NatA catalytic subunit)
Other Name:	
Human Swiss Prot No.:	Q9BSU3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ARD1B at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year