

Antibodies Product Datasheet

Product Name:	ART3 Antibody
Catalog Number:	C49490
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ART3
Concentration:	1 mg/ml
Gene Name:	ART3 TMART
Protein Name:	Ecto-ADP-ribosyltransferase 3 (EC 2.4.2.31) (ADP-ribosyltransferase C2 and C3 toxin-like 3) (ARTC3) (Mono(ADP-ribosyl)transferase 3) (NAD(P)(+)-arginine ADP-ribosyltransferase 3)
Other Name:	
Human Swiss Prot No.:	Q13508
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ART3 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year