

Antibodies Product Datasheet

Product Name:	DHTKD1 Antibody
Catalog Number:	C49494
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DHTKD1
Concentration:	1 mg/ml
Gene Name:	DHTKD1 KIAA1630
Protein Name:	Probable 2-oxoglutarate dehydrogenase E1 component DHKTD1, mitochondrial (EC 1.2.4.2) (Dehydrogenase E1 and transketolase domain-containing protein 1)
Other Name:	
Human Swiss Prot No.:	Q96HY7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DHTKD1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year