

Antibodies Product Datasheet

Product Name:	AGXT Antibody
Catalog Number:	C49526
Reactivity:	Human, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human AGXT
Concentration:	1 mg/ml
Gene Name:	AGXT AGT1 SPAT
Protein Name:	Serine--pyruvate aminotransferase (SPT) (EC 2.6.1.51) (Alanine--glyoxylate aminotransferase) (AGT) (EC 2.6.1.44)
Other Name:	
Human Swiss Prot No.:	P21549
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of AGXT at Human,Rat
Storage and Stability:	Store at -20° C/1 year