

Antibodies Product Datasheet

Product Name:	ENPP4 Antibody
Catalog Number:	C49531
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ENPP4
Concentration:	1 mg/ml
Gene Name:	ENPP4 KIAA0879 NPP4
Protein Name:	Ectonucleotide pyrophosphatase/phosphodiesterase family member 4 (E-NPP 4) (NPP-4) (EC 3.1.-.-)
Other Name:	
Human Swiss Prot No.:	Q9Y6X5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ENPP4 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year