

Antibodies Product Datasheet

Product Name:	HAUS3 Antibody
Catalog Number:	C49547
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HAUS3
Concentration:	1 mg/ml
Gene Name:	HAUS3 C4orf15
Protein Name:	HAUS augmin-like complex subunit 3
Other Name:	
Human Swiss Prot No.:	Q68CZ6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HAUS3 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year