

Antibodies Product Datasheet

Product Name:	CCDC88C Antibody
Catalog Number:	C49583
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CCDC88C
Concentration:	1 mg/ml
Gene Name:	CCDC88C DAPLE KIAA1509
Protein Name:	Protein Daple (Coiled-coil domain-containing protein 88C) (Dvl-associating protein with a high frequency of leucine residues) (hDaple) (Hook-related protein 2) (HkRP2)
Other Name:	
Human Swiss Prot No.:	Q9P219
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CCDC88C at Human, Mouse
Storage and Stability:	Store at -20° C/1 year