

Antibodies Product Datasheet

Product Name:	AP3M1 Antibody
Catalog Number:	C49593
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human AP3M1
Concentration:	1 mg/ml
Gene Name:	AP3M1
Protein Name:	AP-3 complex subunit mu-1 (AP-3 adapter complex mu3A subunit) (Adapter-related protein complex 3 mu-1 subunit) (Mu-adaptin 3A) (Mu3A-adaptin)
Other Name:	
Human Swiss Prot No.:	Q9Y2T2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of AP3M1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year