

Antibodies Product Datasheet

Product Name:	BPGM Antibody
Catalog Number:	C49608
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human BPGM
Concentration:	1 mg/ml
Gene Name:	BPGM
Protein Name:	Bisphosphoglycerate mutase (BPGM) (EC 5.4.2.4) (2,3-bisphosphoglycerate mutase, erythrocyte) (2,3-bisphosphoglycerate synthase) (EC 3.1.3.13) (EC 5.4.2.1) (2,3-diphosphoglycerate mutase) (DPGM) (BPG-dependent PGAM)
Other Name:	
Human Swiss Prot No.:	P07738
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BPGM at Human, Mouse
Storage and Stability:	Store at -20° C/1 year