

Antibodies Product Datasheet

Product Name:	IMP5 Antibody
Catalog Number:	C49613
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human IMP5
Concentration:	1 mg/ml
Gene Name:	SPPL2C IMP5
Protein Name:	Signal peptide peptidase-like 2C (SPP-like 2C) (SPPL2c) (EC 3.4.23.-) (Intramembrane protease 5) (IMP-5)
Other Name:	
Human Swiss Prot No.:	Q8IUH8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IMP5 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year