

Antibodies Product Datasheet

Product Name:	ACY3 Antibody
Catalog Number:	C49624
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ACY3
Concentration:	1 mg/ml
Gene Name:	ACY3 ASPA2
Protein Name:	Aspartoacylase-2 (EC 3.5.1.-) (Acylase III) (Aminoacylase-3) (ACY-3) (Hepatitis C virus core-binding protein 1) (HCBP1) (HCV core-binding protein 1)
Other Name:	
Human Swiss Prot No.:	Q96HD9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ACY3 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year