

Antibodies Product Datasheet

Product Name:	IFT74 Antibody
Catalog Number:	C49647
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human IFT74
Concentration:	1 mg/ml
Gene Name:	IFT74 CCDC2 CMG1
Protein Name:	Intraflagellar transport protein 74 homolog (Capillary morphogenesis gene 1 protein) (CMG-1) (Coiled-coil domain-containing protein 2)
Other Name:	
Human Swiss Prot No.:	Q96LB3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IFT74 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year