

Antibodies Product Datasheet

Product Name:	DDX43 Antibody
Catalog Number:	C49659
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human DDX43
Concentration:	1 mg/ml
Gene Name:	DDX43 HAGE
Protein Name:	Probable ATP-dependent RNA helicase DDX43 (EC 3.6.4.13) (Cancer/testis antigen 13) (CT13) (DEAD box protein 43) (DEAD box protein HAGE) (Helical antigen)
Other Name:	
Human Swiss Prot No.:	Q9NXZ2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of DDX43 at Human
Storage and Stability:	Store at -20° C/1 year