

Antibodies Product Datasheet

Product Name:	EXOSC5 Antibody
Catalog Number:	C49671
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human EXOSC5
Concentration:	1 mg/ml
Gene Name:	EXOSC5 CML28 RRP46
Protein Name:	Exosome complex component RRP46 (Chronic myelogenous leukemia tumor antigen 28) (Exosome component 5) (Ribosomal RNA-processing protein 46) (p12B)
Other Name:	
Human Swiss Prot No.:	Q9NQT4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of EXOSC5 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year