

Antibodies Product Datasheet

Product Name:	IVNS1ABP Antibody
Catalog Number:	C49678
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human IVNS1ABP
Concentration:	1 mg/ml
Gene Name:	IVNS1ABP ARA3 FLARA3 KIAA0850 NS1 NS1BP HSPC068
Protein Name:	Influenza virus NS1A-binding protein (NS1-BP) (NS1-binding protein) (Aryl hydrocarbon receptor-associated protein 3)
Other Name:	
Human Swiss Prot No.:	Q9Y6Y0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IVNS1ABP at Human, Mouse
Storage and Stability:	Store at -20° C/1 year