

Antibodies Product Datasheet

Product Name:	ASCC2 Antibody
Catalog Number:	C49685
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ASCC2
Concentration:	1 mg/ml
Gene Name:	ASCC2 ASC1P100
Protein Name:	Activating signal cointegrator 1 complex subunit 2 (ASC-1 complex subunit p100) (Trip4 complex subunit p100)
Other Name:	
Human Swiss Prot No.:	Q9H1I8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ASCC2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year