

Antibodies Product Datasheet

Product Name:	FAM174A Antibody
Catalog Number:	C49700
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FAM174A
Concentration:	1 mg/ml
Gene Name:	FAM174A NS5ATP6 TMEM157 UNQ1912/PRO4371
Protein Name:	Membrane protein FAM174A (Hepatitis C virus NS5A-transactivated protein 6) (HCV NS5A-transactivated protein 6) (Transmembrane protein 157)
Other Name:	
Human Swiss Prot No.:	Q8TBP5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FAM174A at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year