

Antibodies Product Datasheet

Product Name:	HIVEP1 Antibody
Catalog Number:	C49719
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HIVEP1
Concentration:	1 mg/ml
Gene Name:	HIVEP1 ZNF40
Protein Name:	Zinc finger protein 40 (Cirhin interaction protein) (CIRIP) (Gate keeper of apoptosis-activating protein) (GAAP) (Human immunodeficiency virus type I enhancer-binding protein 1) (HIV-EP1) (Major histocompatibility complex-binding protein 1) (MBP-1) (Posit
Other Name:	
Human Swiss Prot No.:	P15822
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HIVEP1 at Human
Storage and Stability:	Store at -20° C/1 year