

Antibodies Product Datasheet

Product Name:	BCL11A Antibody
Catalog Number:	C49772
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human BCL11A
Concentration:	1 mg/ml
Gene Name:	BCL11A CTIP1 EVI9 KIAA1809 ZNF856
Protein Name:	B-cell lymphoma/leukemia 11A (BCL-11A) (B-cell CLL/lymphoma 11A) (COUP-TF-interacting protein 1) (Ecotropic viral integration site 9 protein homolog) (EVI-9) (Zinc finger protein 856)
Other Name:	
Human Swiss Prot No.:	Q9H165
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BCL11A at Human, Mouse
Storage and Stability:	Store at -20° C/1 year