

Antibodies Product Datasheet

Product Name:	C1orf9 Antibody
Catalog Number:	C49774
Reactivity:	Human, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C1orf9
Concentration:	1 mg/ml
Gene Name:	SUCO C1orf9 CH1 OPT SLP1
Protein Name:	SUN domain-containing ossification factor (Membrane protein CH1) (Protein osteopotential homolog) (SUN-like protein 1)
Other Name:	
Human Swiss Prot No.:	Q9UBS9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of C1orf9 at Human,Rat
Storage and Stability:	Store at -20° C/1 year