

Antibodies Product Datasheet

Product Name:	FRMPD4 Antibody
Catalog Number:	C49784
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FRMPD4
Concentration:	1 mg/ml
Gene Name:	FRMPD4 KIAA0316 PDZD10 PDZK10
Protein Name:	FERM and PDZ domain-containing protein 4 (PDZ domain-containing protein 10) (PSD-95-interacting regulator of spine morphogenesis) (Preso)
Other Name:	
Human Swiss Prot No.:	Q14CM0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FRMPD4 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year