

Antibodies Product Datasheet

Product Name:	IMPG1 Antibody
Catalog Number:	C49786
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human IMPG1
Concentration:	1 mg/ml
Gene Name:	IMPG1 IPM150 SPACR
Protein Name:	Interphotoreceptor matrix proteoglycan 1 (Interphotoreceptor matrix proteoglycan of 150 kDa) (IPM-150) (Sialoprotein associated with cones and rods)
Other Name:	
Human Swiss Prot No.:	Q17R60
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IMPG1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year