

Antibodies Product Datasheet

Product Name:	LRMP Antibody
Catalog Number:	C49787
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human LRMP
Concentration:	1 mg/ml
Gene Name:	LRMP JAW1
Protein Name:	Lymphoid-restricted membrane protein (Protein Jaw1) [Cleaved into: Processed lymphoid-restricted membrane protein]
Other Name:	
Human Swiss Prot No.:	Q12912
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LRMP at Human
Storage and Stability:	Store at -20° C/1 year