

Antibodies Product Datasheet

Product Name:	KANSL2 Antibody
Catalog Number:	C49848
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KANSL2
Concentration:	1 mg/ml
Gene Name:	KANSL2 C12orf41 NSL2
Protein Name:	KAT8 regulatory NSL complex subunit 2 (NSL complex protein NSL2) (Non-specific lethal 2 homolog)
Other Name:	
Human Swiss Prot No.:	Q9H9L4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KANSL2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year