

Antibodies Product Datasheet

Product Name:	LDB1 Antibody
Catalog Number:	C49850
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human LDB1
Concentration:	1 mg/ml
Gene Name:	LDB1 CLIM2
Protein Name:	LIM domain-binding protein 1 (LDB-1) (Carboxyl-terminal LIM domain-binding protein 2) (CLIM-2) (LIM domain-binding factor CLIM2) (hLdb1) (Nuclear LIM interactor)
Other Name:	
Human Swiss Prot No.:	Q86U70
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LDB1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year