

Antibodies Product Datasheet

Product Name:	KIR2DS4 Antibody
Catalog Number:	C49866
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KIR2DS4
Concentration:	1 mg/ml
Gene Name:	KIR2DS4 CD158I KKA3 NKAT8
Protein Name:	Killer cell immunoglobulin-like receptor 2DS4 (CD158 antigen-like family member I) (MHC class I NK cell receptor) (Natural killer-associated transcript 8) (NKAT-8) (P58 natural killer cell receptor clones CL-39/CL-17) (p58 NK receptor CL-39/CL-17) (CD ant
Other Name:	
Human Swiss Prot No.:	P43632
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KIR2DS4 at Human
Storage and Stability:	Store at -20° C/1 year