

Antibodies Product Datasheet

Product Name:	ENPP6 Antibody
Catalog Number:	C49899
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ENPP6
Concentration:	1 mg/ml
Gene Name:	ENPP6 UNQ1889/PRO4334
Protein Name:	Ectonucleotide pyrophosphatase/phosphodiesterase family member 6 (E-NPP 6) (NPP-6) (EC 3.1.4.-) (EC 3.1.4.38) (Choline-specific glycerophosphodiester phosphodiesterase) (Glycerophosphocholine cholinephosphodiesterase) [Cleaved into: Ectonucleotide pyropho
Other Name:	
Human Swiss Prot No.:	Q6UWR7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ENPP6 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year