

Antibodies Product Datasheet

Product Name:	GPKOW Antibody
Catalog Number:	C49909
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GPKOW
Concentration:	1 mg/ml
Gene Name:	GPKOW GPATC5 GPATCH5 T54
Protein Name:	G patch domain and KOW motifs-containing protein (G patch domain-containing protein 5) (Protein MOS2 homolog) (Protein T54)
Other Name:	
Human Swiss Prot No.:	Q92917
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GPKOW at Human, Mouse
Storage and Stability:	Store at -20° C/1 year