

Antibodies Product Datasheet

Product Name:	CPVL Antibody
Catalog Number:	C49913
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CPVL
Concentration:	1 mg/ml
Gene Name:	CPVL VLP PSEC0124 UNQ197/PRO223
Protein Name:	Probable serine carboxypeptidase CPVL (EC 3.4.16.-) (Carboxypeptidase, vitellogenic-like) (Vitellogenic carboxypeptidase-like protein) (VCP-like protein) (hVLP)
Other Name:	
Human Swiss Prot No.:	Q9H3G5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CPVL at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year