

Antibodies Product Datasheet

Product Name:	HMGXB4 Antibody
Catalog Number:	C49917
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HMGXB4
Concentration:	1 mg/ml
Gene Name:	HMGXB4 HMG2L1 HMGB CG
Protein Name:	HMG domain-containing protein 4 (HMG box-containing protein 4) (High mobility group protein 2-like 1) (Protein HMGB CG)
Other Name:	
Human Swiss Prot No.:	Q9UGU5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HMGXB4 at Human
Storage and Stability:	Store at -20° C/1 year