

Antibodies Product Datasheet

Product Name:	ASAH2 Antibody
Catalog Number:	C49941
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ASAH2
Concentration:	1 mg/ml
Gene Name:	ASAH2 HNAC1
Protein Name:	Neutral ceramidase (N-CDase) (NCDase) (EC 3.5.1.23) (Acylsphingosine deacylase 2) (BCDase) (LCDase) (hCD) (N-acylsphingosine amidohydrolase 2) (Non-lysosomal ceramidase) [Cleaved into: Neutral ceramidase soluble form]
Other Name:	
Human Swiss Prot No.:	Q9NR71
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ASAH2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year