

Antibodies Product Datasheet

Product Name:	EXOSC9 Antibody
Catalog Number:	C49956
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human EXOSC9
Concentration:	1 mg/ml
Gene Name:	EXOSC9 PMSCL1
Protein Name:	Exosome complex component RRP45 (Autoantigen PM/Scl 1) (Exosome component 9) (P75 polymyositis-scleroderma overlap syndrome-associated autoantigen) (Polymyositis/scleroderma autoantigen 1) (Polymyositis/scleroderma autoantigen 75 kDa) (PM/Scl-75)
Other Name:	
Human Swiss Prot No.:	Q06265
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of EXOSC9 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year