

Antibodies Product Datasheet

Product Name:	ARH3 Antibody
Catalog Number:	C49957
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ARH3
Concentration:	1 mg/ml
Gene Name:	ADPRHL2 ARH3
Protein Name:	Poly(ADP-ribose) glycohydrolase ARH3 (EC 3.2.1.143) (ADP-ribosylhydrolase 3) ([Protein ADP-ribosylarginine] hydrolase-like protein 2)
Other Name:	
Human Swiss Prot No.:	Q9NX46
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ARH3 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year