

Antibodies Product Datasheet

Product Name:	C17orf80 Antibody
Catalog Number:	C49966
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C17orf80
Concentration:	1 mg/ml
Gene Name:	C17orf80 HLC8 MIG3
Protein Name:	Uncharacterized protein C17orf80 (Cell migration-inducing gene 3 protein) (Human lung cancer oncogene 8 protein) (HLC-8)
Other Name:	
Human Swiss Prot No.:	Q9BSJ5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of C17orf80 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year