

Antibodies Product Datasheet

Product Name:	FCRLA Antibody
Catalog Number:	C50008
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FCRLA
Concentration:	1 mg/ml
Gene Name:	FCRLA FCRL FCRL1 FCRLM1 FCRX FREB UNQ291/PRO329
Protein Name:	Fc receptor-like A (Fc receptor homolog expressed in B-cells) (Fc receptor-like and mucin-like protein 1) (Fc receptor-like protein) (Fc receptor-related protein X) (FcRX)
Other Name:	
Human Swiss Prot No.:	Q7L513
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FCRLA at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year