

Antibodies Product Datasheet

Product Name:	C20orf4 Antibody
Catalog Number:	C50022
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C20orf4
Concentration:	1 mg/ml
Gene Name:	AAR2 C20orf4 CGI-23 PRO0225
Protein Name:	Protein AAR2 homolog (AAR2 splicing factor homolog)
Other Name:	
Human Swiss Prot No.:	Q9Y312
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of C20orf4 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year