

Antibodies Product Datasheet

Product Name:	HLA-DQA2 Antibody
Catalog Number:	C50034
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human HLA-DQA2
Concentration:	1 mg/ml
Gene Name:	HLA-DQA2 HLA-DXA
Protein Name:	HLA class II histocompatibility antigen, DQ alpha 2 chain (DX alpha chain) (HLA class II histocompatibility antigen, DQ(6) alpha chain) (HLA-DQA1) (MHC class II DQA2)
Other Name:	
Human Swiss Prot No.:	P01906
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of HLA-DQA2 at Human
Storage and Stability:	Store at -20° C/1 year