

Antibodies Product Datasheet

Product Name:	IMPAD1 Antibody
Catalog Number:	C50035
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human IMPAD1
Concentration:	1 mg/ml
Gene Name:	IMPAD1 IMPA3
Protein Name:	Inositol monophosphatase 3 (IMP 3) (IMPase 3) (EC 3.1.3.25) (EC 3.1.3.7) (Golgi 3-prime phosphoadenosine 5-prime phosphate 3-prime phosphatase) (Golgi-resident PAP phosphatase) (gPAPP) (Inositol monophosphatase domain-containing protein 1) (Inositol-1(or
Other Name:	
Human Swiss Prot No.:	Q9NX62
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IMPAD1 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year