

Antibodies Product Datasheet

Product Name:	LRRC29 Antibody
Catalog Number:	C50040
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human LRRC29
Concentration:	1 mg/ml
Gene Name:	LRRC29 FBL9 FBXL9
Protein Name:	Leucine-rich repeat-containing protein 29 (F-box and leucine-rich repeat protein 9) (F-box protein FBL9) (F-box/LRR-repeat protein 9)
Other Name:	
Human Swiss Prot No.:	Q8WV35
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LRRC29 at Human
Storage and Stability:	Store at -20° C/1 year