

Antibodies Product Datasheet

Product Name:	ARHGEF18 Antibody
Catalog Number:	C50092
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ARHGEF18
Concentration:	1 mg/ml
Gene Name:	ARHGEF18 KIAA0521
Protein Name:	Rho guanine nucleotide exchange factor 18 (114 kDa Rho-specific guanine nucleotide exchange factor) (p114-Rho-GEF) (p114RhoGEF) (Septin-associated RhoGEF) (SA-RhoGEF)
Other Name:	
Human Swiss Prot No.:	Q6ZSZ5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ARHGEF18 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year